Updated 5/24

XXXXXXXXXX 40ml bottles

) box of 24 (

**WE NO LONGER Supply

Needles/Holders/Phlebotomy Supplies**

Please reach out to warehouse or your

office supply vendor to order

() box of 24 () case of 96

) case of 96

Please fill-in quantity to the left box of the item needed

MML Supply Request Form

Phone: 301-665-4900

Custom

Custom

Surgical Path Reqs (blue)

Non-Gyn cytology forms

Courier labels (Medspeed)

Patient labels (30 per sheet)

Epic printer labels

ABN forms

1 part padded

2 part duplicate

Fax: 301-665-4942

Urine for Cytlology(Urology CTR)

Dark blue lidded cups should **NOT

be used for urine specimens**

thin-prepcytolyt solution

Special Request Items (not listed above)

Hereditary Cancer kit (Quest)

Date: Phone:					Orders are filled on Thursday mornings				
Practice Name: Contact person:				Deadline		Deadline to	to fax orders is Wednesday evening		
Address:					Fax to		301-665-4942		
Please in	dicate the Quantity to have	your order fi	illed						
Quantity	PHLEBOTOMY SUPPLIES	Quantity	SWABS	Quantity	GENERAL S	UPPLIES	Quantity	PAP SUPPLIES	
	Gold tubes (SST)		White (DRS/SSCR/WOUND)		Transport bag	S		Thin prep vials	
	Lavender tubes		Dk Green (eye/ear)		Small			tray of 25	
	Light blue tubes (Na Citrate)		Royal Blue (Bordetella)		Large Stat			case of 250	
	Green top (Na Hep)		Light Blue top (MRSA/Herpes)					Brushes/spatulas	
	Green top (Lithium Hep)		BASC (Tric/Candida/Gardnerella)					Brooms	
	Light Green top (PST & Gel Lith.Hep)		Flu/RSV/Covid swabs UTM 3ml					Pap forms	
	Microtainer tubes (specify below)		Other Viruses Culture/UTM 3ml						
	lavender	XXXXXXXXX	стдс				Quantity STOOL COLLECTION		
	amber-SST		Urine (male/female)	Quantity	URINE COL	LECTION		Culture (green top)	
	Dark blue (Serum)		Swab (female only)		Vacutainer Cu	ps (with needle)		Ova & parasite	
	Dark blue (EDTA)	Quantity	FORMS & LABELS		Vacutainer Kits (cup & tubes)			Occult blood cards	
	Pink tubes (Blood bank)	XXXXXXXXX	Lab Slips-please send copy with order and		Gray tubes (ur	rine culture)		Sterile cups	
	Red top tubes		make sure quanity is indicated		24 hour collec	tion jug (regular)			
XXXXXXXXX	Blood Culture bottles	Questions or Issues, Please contact: Melissa Benefield 301-665-4929 OR Rose Monn 301-665-4932			24 hr trace ele	ement jug (for metals)	Quantity	MISCELLANEOUS	
	adult				Towelettes (W	/ipes)		Fertility bottles	
	pediatric		Generic 1 part padded		Pediatric U-bags			Sputum collection cups	
Quantity	BIOPSY SUPPLIES		Generic 2 part duplicate		Hats			Autoclave kits (KILITS)	